

# GABA rabbit pAb – IS1006

Ref: IS1006-sp

The anti-GABA antibody IS1006 rabbit polyclonal antibody displays high affinity and specificity and allows direct detection of GABA in whole mounts, cell culture and tissue sections when samples are prepared using the [STAINperfect immunostaining kit A](#).

<b>Clonality</b>	Polyclonal antibody
<b>Host</b>	Rabbit
<b>Reactivity</b>	Reacts with all species
<b>Tested samples</b>	Whole mounts, cell culture, tissue sections
<b>Staining procedure</b>	<a href="#">STAINperfect immunostaining kit A</a>
<b>Format</b>	50µL (approx. 40 tissue sections)
<b>References</b>	<a href="#">Citations in literature</a>

## INFORMATIONS

---

### Product overview

---

<b>Product name</b>	GABA antibody - Rabbit Polyclonal antibody
<b>Synonyms</b>	Anti-Gamma-Aminobutyric acid antibody
<b>Immunogen</b>	Conjugated GABA
<b>Specificity</b>	When tested in competitive ELISA, the anti-conjugated GABA antibody did not show any significant cross reactivity with Gamma-Aminobutyric acid analogs, including Beta-Alanine and D-Alanine
<b>Volume</b>	50µL

### Storage

---

<b>Form</b>	Liquid
<b>Purity</b>	Purified anti-serum
<b>Storage</b>	Store at +4°C for short term (1-2 months). Aliquot and store at -20°C for long term. Avoid repeated freeze / thaw cycles
<b>Material safety datasheet</b>	<a href="#">Download MSDS</a>

# PROTOCOLS


**IF - Cell cultures, Whole mounts, Tissue sections** Dilute antibody with the antibody diluent provided in the [STAINperfect immunostaining kit A](#). Use at 1/250 -1/1000 dilution. Follow the STAINperfect protocol suited to your sample

**Comments** Optimal working dilutions must be determined by the end-user

**Restrictions** For research use only

**Full protocol** [Download STAINperfect protocol for GABA staining](#)

## Protocols-at-a-glance



[Complete Instructions for Use](#)

[Protocol-at-a-glance for cell cultures](#)

[Protocol-at-a-glance for whole mounts](#)

[Protocol-at-a-glance for tissue sections](#)

## REFERENCES

---

### Product citations:


- [Villamayor et al. Structural, morphometric and immunohistochemical study of the rabbit accessory olfactory bulb. Brain Struct Funct. 2019 Dec 4.](#)
- [Smits et al., Single-cell transcriptomics reveals multiple neuronal cell types in human midbrain-specific organoids. bioRxiv. 2019, March 28.](#)
- [Yokoi et al. Impact of Sleep-Wake-Associated Neuromodulators and Repetitive Low-Frequency Stimulation on Human iPSC-Derived Neurons Front Neurosci. 2019 May 29.](#)
- [Traub et al. hiPS Cell-Derived Neurons for High-Throughput Screening. Methods Mol Biol. 2019;1994:243-263](#)

### Selected articles on GABA:

- [Liu X, Wang Q, Haydar TF, Bordey A. Nonsynaptic GABA signaling in postnatal subventricular zone controls proliferation of GFAP-expressing progenitors. Nat Neurosci. 2005 Sep;8\(9\):1179-87. Epub 2005 Aug 14.](#)
- [Lawrence JJ. Cholinergic control of GABA release: emerging parallels between neocortex and hippocampus. Trends Neurosci. 2008 Jul;31\(7\):317-27. doi: 10.1016/j.tins.2008.03.008. Epub 2008 Jun 13.](#)
- [Baulac S, Huberfeld G, Gourfinkel-An I, Mitropoulou G, Beranger A, Prud'homme JF, Baulac M, Brice A, Bruzzone R, LeGuern E. First genetic evidence of GABA\(A\) receptor dysfunction in epilepsy: a mutation in the gamma2-subunit gene. Nat Genet. 2001 May;28\(1\):46-8.](#)

## Product pictures


---



### GABA labeling (green) in mouse cortical primary neurons


Immunodetection of GABA- (green) and MAP2- (red) positive neurons in mouse primary cortical culture. GABAergic neurons were stained using IS1006 anti-GABA antibody combined with MAP2 antibody using the stainperfect immunostaining kit A and appropriate protocol, respectively. Alexa fluor 488 and 546 secondary antibodies were used and pictures were captured by high-content imaging.

---




### Affinity & specificity of anti-GABA antibody

Competitive ELISA demonstrates that low amounts of GABA conjugate are required to abolish antigen-antibody reaction (high affinity), while rising concentrations of  $\beta$ -Alanine and D-Alanine conjugates do not affect reaction (high specificity).



### L-Glutamate and GABA in adult mouse primary cortical neurons


Adult mouse primary cortical neurons were stained with mouse monoclonal anti-L-glutamate antibody (red) combined with anti-GABA rabbit polyclonal antibody (green). Staining was performed using optimized sample preparation with STAINperfect immunostaining kit A and according to the protocol for cell culture. Fluorescent labeled secondary antibody were used and pictures were acquired by confocal imaging.



### GABA labeling (green) in rat cortical primary neurons

Immunodetection of GABA- (green) and MAP2- (red) positive neurons in rat primary cortical culture. Total and GABAergic neurons were stained using an anti-MAP2 antibody and IS anti-GABA antibody combined with the stainperfect immunostaining kit A and appropriate protocol, respectively. Alexa fluor 546 and 488 secondary antibodies were used and pictures were captured by high-content imaging.

### Gamma-aminobutyric acid (GABA)



In the mammalian brain, inhibitory neurotransmitter Gamma-aminobutyric acid (GABA) is mainly synthesized from excitatory L-Glutamate by enzyme glutamic acid decarboxylase (GAD). Regulating neuronal excitability, GABAergic synapses are present throughout the CNS, although GABA is found most highly concentrated in the substantia nigra, the globus pallidus nuclei, the hypothalamus, the periaqueductal grey matter and the hippocampus.

## Contact information

Immusmol  
229 Cours de l'Argonne  
33 000 Bordeaux - France  
Tel: +33 (0) 5 6431 1170  
www.immusmol.com

**To order, review, ask for technical support, visit product page at:**

<https://www.immusmol.com/shop/gaba-polyclonal-antibody-bundle/>