

# Androstenediol-Glucuronide ELISA kit

## I Fast I Serum

Ref: AA-E-1500

This Androstenediol-Glucuronide ELISA kit enables the quantitative measurement of C19 steroid 3 $\alpha$ -Diol-G in serum samples to evaluate androgen metabolism. Intended for *in vitro* diagnostics and research applications, the test comes in 96-well plate format, allowing the screening of large experimental series within short lead times.

<b>Sample type</b>	Serum
<b>Capacity</b>	96 tests
<b>Sensitivity</b>	0.1 ng/mL
<b>Range</b>	0/0.25 - 50 ng/mL
<b>Assay time<sup>75</sup></b>	45 min
<b>Reactivity</b>	Reacts with all species

## INFORMATIONS

---

<b>Product name</b>	Androstenediol-Glucuronide ELISA
<b>Description</b>	<i>In vitro</i> diagnostics / Research enzyme immunoassay (ELISA) for the quantitative determination of Androstenediol-Glucuronide ELISA in serum samples
<b>Labels</b>	IVD (EU only) or RUO, CE
<b>Format</b>	96-well plate
<b>Samples</b>	Serum
<b>Sample volume</b>	50µL
<b>Reactivity</b>	Reacts with all species
<b>Standard range</b>	0/0.25 – 50 ng/mL
<b>Sensitivity</b>	0.1 ng/mL
<b>Specificity</b>	No significant cross-reactivity was observed with 3α-Diol-G analogs such as Testosterone, Progesterone, Androstenedione, and Cortisol.
<b>Assay time</b>	45 min

## PROTOCOLS

---

<b>Product storage</b>	Store at 2-8°C for up to 6 months
<b>Sample collection &amp; storage</b>	<u>Serum</u> : Store samples at 2-8°C for up to 48h or -20°C for longer period (up to 6 months).
<b>Sample preparation</b>	Enzymatic conjugation (30 min)
<b>ELISA</b>	Incubation, revelation and read steps (15min)
<b>Detailed protocol</b>	IVD: <a href="#">Download instructions for use</a> RUO: <a href="#">Download instructions for use</a>

## Product pictures

---



### Androstanediol-Glucuronide ELISA kit AA-E-1500

Androstanediol-Glucuronide ELISA kit AA-E-1500

---

## Contact information

---

Immusmol  
229 Cours de l'Argonne  
33 000 Bordeaux - France  
Tel: +33 (0) 5 6431 1170  
[www.immusmol.com](http://www.immusmol.com)

---

**To order, review, ask for technical support, visit product page at:**

<https://www.immusmol.com/shop/androstanediol-glucuronide-elisa-kit-i-fast-i-serum/>